

coretech™

vive health innovation

Coretech 916 Wrist Brace

SUP1039

HCPC: L3915 - L3916

The 916 Wrist Brace provides compression and support during the healing process of common wrist injuries. Crafted with a light and breathable materials and a bendable inner metal splint that conforms to the shape of your wrist. This brace has a universal fit and can be used for either the left or right hand.

Specifications

- What's Included:
 - 1pc 916 Wrist Brace
- Universal Size, fits left and right adult arm
- Color:
 - Black
- Materials:
 - Terylene
 - Aluminum
 - Nylon
 - PVC
 - Latex-free
- Dimensions:
 - 11.6 in (29.5 cm) Length
- Weight:
 - 0.53 lb (0.24 kg)

Distributed by

**vive
health**

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com





Premium Brace Quality

Created from lightweight perforated material that stays comfortable and resists excess moisture and heat during all-day use.



Customizable Design

Easily adjust the ROM hinge for a personalized level of compression and flexion that can be changed according to your needs.



Promotes Healing

Provides increased circulation and wrist support that may help the recovery and relief of carpal tunnel syndrome, osteoarthritis, rheumatoid arthritis, and other common wrist conditions.



Universal Sizing

Ambidextrous fit suits men and women of all sizes. The flexible brace naturally conforms to a variety of wrist and forearm shapes.

Distributed by

**vive
health**

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

☎ (239) 220-5772 🌐 vivehealth.com ✉ service@vivehealth.com

📍 8955 Fontana Del Sol Way, Naples, FL 34109, USA