

Coretech - 350 Ankle Air Splint

SUP2043

HCPCS: L4350

Providing cushioning support for injured or weak ankles, the Coretech 350 Ankle Air Splint utilizes therapeutic compression to relieve pain and increase circulation while reducing inflammation and swelling in the ankle following surgery or injury. The exterior lateral and medial supports also stabilize the ankle for safe, easy mobility throughout the day.

Universally designed for the right or left ankle, the ankle air splint features a slim heel strap that can be lengthened for a personalized fit. This also allows the brace to be used comfortably with any shoe style or size. A small hand pump is included for easily inflating and deflating the side air panels for customizable compression and optimal ankle support.

Specifications

- What's Included:
 - 1pc 350 ankle air splint
- 1pc air pump
- Colors:
 - Black outer shell and green inner padding
- Sizes:
 - Universal size
- Materials:
 - PP plastic, Nylon, Iron, Polyester, PVC, 100% PU, Rubber, Latex-free
- Care Instructions:
 - Wipe down with damp cloth, dry immediately, the pads can be removed from each side.
- Sizing Information: The bottom is adjustable so it can fit men and women
 - Men 4 - 13+
 - Women 7 - 15+

Distributed by

**vive
health**

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com





Rigid Side Supports

Designed to provide stabilizing support for ankle sprains, injury and instability.



Customizable Compression

Inflatable side panels stabilize the ankle and increase circulation with therapeutic compression.



Easy Inflation

Includes a small hand pump to adjust the pre-inflated side air panels.



Works With Any Shoe Style

Low profile heel strap allows the brace to be worn with any shoe style for all-day support.

Distributed by

**vive
health**

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

(239) 220-5772  vivehealth.com  service@vivehealth.com

 8955 Fontana Del Sol Way, Naples, FL 34109, USA