

# coretech™

vive health innovation

## Coretech ROM Knee Brace

SUP2071

HCPC: L1833/L1832

Promote healing and help prevent further injury by stabilizing and supporting the knee. The Coretech ROM Knee Brace fits most adult leg sizes with an extension range between 18" to 23.5" and adjustable fastening straps. Anti-slip lining and aluminum alloy backing make wearing over or under clothing comfortable and light. Customize the range of motion with the easy-to-use locking dial. A controlled range of motion between 120° to open flexion and 90° to 0° of extension helps provide the support needed for a more comfortable and stable healing process. The flexibility of motion assists with every stage of recovery. Ideal for following strains, patella realignment, ACL repair, and other knee-related surgeries.

### Specifications

- What's Included:
  - 1 ROM Knee Brace
- Universal Size, fits left and right adult leg.
- Color:
  - Black
- Materials:
  - Aluminium alloy,
  - PP plastic,
  - Nylon,
  - Polyester.
  - Latex-free.
- Dimensions:
  - 18.5 x 7 x 3.5 in /
  - 47 x 17.8 x 8.9 cm (L x W x H)
- Weight:
  - 2 lb / 0.9 kg
- Length range:
  - 18 – 23.5 in / 45.7 – 56.7 cm
- Adjustable Flexion:
  - 90°, 70°, 40°, 30°
- Adjustable Extension:
  - 30°, 10°, 0°, -10°

Distributed by

**vive  
health**

8955 Fontana Del Sol Way  
Naples, FL 34109  
1-800-487-3808  
[www.vivehealth.com](http://www.vivehealth.com)





## Superior Knee Recovery Support

Helps reinforce and stabilize the knee to reduce re-injury risks and encourage proper healing.



## Reversible & Accommodating Sizing

Compatible with both left and right legwear. With non-slip fastening straps that can be cut for further customization.



## Lightweight & Secure Fit

Aluminum alloy construction allows a weightless feel. Easy to open and release buckles makes wearing over clothes simple and comfortable.



## Adjustable Range Of Motion

Fully customize the degree of movement. Use the locking dial to set the preferred angle in place. Extension of 30° to -10° and flexion of 30° to 90°.

Distributed by

**vive**  
**health**

8955 Fontana Del Sol Way  
Naples, FL 34109  
1-800-487-3808  
[www.vivehealth.com](http://www.vivehealth.com)

☎ (239) 220-5772    🌐 [vivehealth.com](http://vivehealth.com)    ✉ [service@vivehealth.com](mailto:service@vivehealth.com)

📍 8955 Fontana Del Sol Way, Naples, FL 34109, USA