



vive health innovation

Coretech 762 Elbow Brace

SUP2090

HCPC: L3762

Coretech 762 Elbow Brace holds the elbow joint at a 106° - 110° angle to treat and manage pain from cubital tunnel syndrome, ulnar nerve entrapment, elbow hyperextension, and surgery recovery.



Specifications

- What's Included:
 - 762 Cubital Tunnel Elbow Brace (left or right)
- Color:
 - Black
- Weight:
 - 0.52 lb. / 0.62 kg.
- Materials:
 - ABS
 - Polyurethane
 - Polyester
 - Nylon
- Latex-free:
 - Yes

SIZE CHART

\	Top Max. Circum.	Bottom Max. Circum.
S / M	9.84" to 21.06" (25 to 53.5 cm)	8.66" to 2.99" (22 to 33 cm)
L / XL	10.8" to 21.5" (27.5 to 54.5 cm)	8" to 13.8" (20.5 to 35 cm)

Distributed by



8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com



Pain-relieving Support

Immobilizes and positions the elbow at 106°-110° to help reduce pain and numbness associated with Cubital Tunnel Syndrome.



Versatile Fit

Suitable for both men and women, customize fit on left or right arm.



Lightweight And Breathable

Constructed with a high-density shell and perforated liner with adjustable pads for greater protection and comfort.



Adjustable

3 adjustable straps and easy clip-on buckles to make application as easy as possible.

Distributed by

**vive
health**

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

☎ (239) 220-5772 🌐 vivehealth.com ✉ service@vivehealth.com
📍 8955 Fontana Del Sol Way, Naples, FL 34109, USA