

## Coretech™ 637 Lumbar Brace

SUP2026

HCPCS: L0637 / L0650

Providing the maximum amount of support, the Coretech 637 LSO Lumbar Brace utilizes an elongated, rigid back panel and built-in side panels to stabilize the trunk, complementing current therapies and promoting healing following injury or surgery. A removable pulley system allows the level of targeted compression to be easily adjusted for optimal support and patient comfort while participating in everyday activities.

The brace effectively relieves pain and pressure due to muscle fatigue, strains, scoliosis, herniated discs, sciatica and spondylolisthesis. Lightly padded, the brace also includes an adjustable waistband, making it easier to fit to the individual, promoting higher patient compliance rates.

### Specifications

- What's Included
  - 1x 637 Lumbar Brace
- Colors: Black
- Adjustable size: S-2XL (25" - 50" waist)
- Dimensions: 15" x 9.5" x 5.5"
  - Back splint: 15"(H) x 10" (W)
  - Side splints: 7" x 7.5"
- 1.76 lbs
- Latex Free

Materials:  
Nylon oxford cloth 10%, brush cloth 10%, foam 20%, ABS plastic 40%, PP plastic 20%

Care:  
Care: Spot clean with a damp cloth, air dry. Do NOT machine wash or dry.

Distributed by

**vive  
health**

8955 Fontana Del Sol Way  
Naples, FL 34109  
1-800-487-3808  
[www.vivehealth.com](http://www.vivehealth.com)





### Rigid Splints

Includes a large rigid back panel and side splints for exceptional lumbar support and trunk stability.

### Adjustable Waist Straps

Easily adjust the side panels of the brace to fit waist sizes ranging from 25" to 50".

### Comfort Blend

Lightly padded for greater comfort, the brace is made with a soft, lightweight fabric blend.

### Removable Pulley System

Featuring wide pull tabs, the brace utilizes a secure pulley system for customizable compression and support.



### Premium Fastener System

Easily secured with strong, dual fasteners for a personalized, non-slip fit.



Distributed by

**vive health**  
 8955 Fontana Del Sol Way  
 Naples, FL 34109  
 1-800-487-3808  
 www.vivehealth.com

☎ (239) 220-5772    🌐 vivehealth.com    ✉ service@vivehealth.com  
 📍 8955 Fontana Del Sol Way, Naples, FL 34109, USA