



vive health innovation

CORETECH 807 Thumb Splint

SUP2040/SUP2041

HCPCS: L3807/L3809

Immobilizing the thumb in a neutral position, the Coretech 807 thumb splint provides pain-relieving support for sprains, strains and scaphoid injuries. The spica splint includes a rigid aluminum splint that is removable for greater flexibility and two malleable splints for customizable support.

Lightweight and breathable, the thumb brace is constructed with a soft polyester blend and a memory foam lining for exceptional comfort with four adjustable fastener straps to provide therapeutic compression and a secure fit. The splint also allows the fingers to move freely to easily complete daily tasks, encouraging a better rate of patient compliance.



Specifications

- ✦ What's Included:
 - 1 pc. Wrist support (L or R)
- ✦ Colors:
 - Black with blue interior
- ✦ Sizes:
 - XS - XL
- ✦ Materials:
 - Nylon
 - PP
 - Metal
 - Aluminum
 - Sponge
 - PU
 - Velvet
 - Polyester
 - Cotton
 - Latex free
- ✦ Care Instructions:
 - Spot clean, do not machine wash or dry. Keep dry.
- ✦ Sizing Information: XS 4.5 - 5.5 in.
 - S 5.5 - 6.6 in.
 - M 6.5 - 7.5 in.
 - L 7.5 - 8.5 in.
 - XL Over 8.5 in.

Distributed by



8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com



Immobilizing Support

Limiting movement of the thumb, the thumb brace includes one rigid splint that is removable when more flexibility is needed.



Flexible Side Splints

Providing customizable support, the thumb brace includes two malleable side splints.



Breathable Design

Constructed with a lightweight fabric blend, the padded thumb brace features ventilation holes for breathability and comfort.



Adjustable Fit

Each strap is easily adjustable for a personalized fit.

Distributed by

**vive
health**

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

☎ (239) 220-5772 🌐 vivehealth.com ✉ service@vivehealth.com

📍 8955 Fontana Del Sol Way, Naples, FL 34109, USA