

coretech™
vive health innovation

Coretech™ 396 Dorsal Night Splint

SUP2045

HCPCS: L4396 / L4397

Perfect for side, back and stomach sleepers, the Coretech 396 Dorsal Night Splint gently stretches the plantar fascia and soft tissues of the foot and calf while supporting proper foot positioning to effectively relieve pain and reduce inflammation. The rigid composite night splint is ideal for relieving plantar fasciitis, foot drop, achilles tendonitis, heel pain and calf muscle cramps.

The open design is lightweight and breathable for greater comfort throughout the night and is completely padded to eliminate chafing and hot spots on the foot, ankle or calf. Easy to adjust, the dorsal night splint can be worn on the right or left foot and includes a heel strap to ensure a secure, non slip fit even on the most restless nights.

Specifications

- What's Included
 - 1x 396 Dorsal Night Splint
- Color: Black
- Sizes: S/M, L/XL
- Dimensions: 13" x 6" x 4.5"
- Weight: 0.54lbs
- Latex Free

Materials:
Nylon, polyester, ABS, neoprene

Care:
Remove plastic splint, hand wash with cold water and mild soap. Air Dry. Do NOT machine wash or dry.



Distributed by

**vive
health**

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com



Low Profile Design

Featuring an open toe and heel, the low profile splint is comfortable for side, stomach or back sleepers.

Comfortable Cushioning

Each section of the adjustable night splint is lined with additional padding to eliminate chafing and irritation.

Lightweight & Breathable

Made with ultra-soft polyester and neoprene that is lightweight and breathable for exceptional nighttime comfort.

Rigid Nighttime Support

Gently stretches and properly positions the foot and ankle throughout the night to relieve pain due to plantar fasciitis, soft tissue damage, arthritis and foot drop..



Premium Fastener System

Easily secure with strong fasteners for a personalized, non slip fit.

SIZE	HEIGHT	MAX CALF CIRCUMFERENCE	WEIGHT
S/M	9"	15"	7 oz.
L/XL	11"	19"	10 oz.

Distributed by



8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

☎ (239) 220-5772 🌐 vivehealth.com ✉ service@vivehealth.com

📍 8955 Fontana Del Sol Way, Naples, FL 34109, USA