

coretech™

vive health innovation

Coretech 833 Knee Brace

SUP2062

HCPCS: L1833 / L1832

The Coretech 833 Knee Brace offers support for knee and meniscus injuries such as MCL and ACL tears (medial collateral ligament and anterior cruciate ligament), strains and sprains, while also reducing the risk of re-injury. Providing medical and lateral support, a hinge angle immobilization selector, above and below knee support straps, and mid-knee compression, the Coretech 833 Knee Brace is the ultimate injury recovery brace. It is recommended to adjust the hinge angle before putting the brace on, however, hinge adjustments may be made while wearing the brace if preferred.



Specifications

- What's Included:
 - Adjustable hinged knee support with straps and buckles
- Colors:
 - Black
- Sizes:
 - Small, Medium, Large
- Materials:
 - Neoprene, Nylon, Polyester, Spandex, foam
- Care Instructions:
 - Do not machine wash or dry. Spot clean with a damp cloth and air dry.
- Sizing Information:
 - Small: Up to 19"
 - Medium: 19" - 23"
 - Large: 23" - 29"

Distributed by

**vive
health**

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com



Anti-Slip Bands

Designed with internal silicone bands to prevent slipping and shifting during normal and heavy-duty activity. It is easy to put on and lightweight for non-contact sports.



Double Hinged

Dual hinge design with extension and flexion aluminum pins for easy on-the-go adjustments. Allowing the knee to bend from 0, 10, 20, 40, 60, 90 or 150 degrees.



Highly Adjustable

Integrated easy clip-on buckles and trimmable straps for the perfect fit. It helps progress through injury recovery or post-surgery rehab.



Breathable

High-performance neoprene promotes air circulation throughout the knee and it is latex-free.

Distributed by

**vive
health**

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

☎ (239) 220-5772 🌐 vivehealth.com ✉ service@vivehealth.com

📍 8955 Fontana Del Sol Way, Naples, FL 34109, USA