

Coretech 621 SI Belt

SUP2057BLK

HCPCS: L0621 / L0622

Comfortable for all-day wear, the Coretech 621 SI Belt is great to use for maximum support at home, at the office, and for play.

This belt compresses and supports sacroiliac joints, helping to treat muscle pain and soreness.

Specifications

- What's Included:
 - S/M SI belt
 - L/XL SI belt
 - 8" Velcro extender
- Colors:
 - Black
- Sizes:
 - S/M - L/XL
- Materials:
 - 50% polyester and 50% nylon
 - Latex free
- Care Instructions:
 - Spot clean with a damp cloth, do NOT machine wash or dry
- Sizing Information:
 - S/M: Up to 32" waist circumference
 - L/XL: 32" - 52"
 - 8" extender: adds 8" to any belt size

Distributed by

**vive
health**

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com





Adjustable And Secure Fit

Helps stabilize and reinforce the SI joint during daily activities. Attached pulley system on the back allows user to adjust the compression levels as needed and provides rigid back support. Thumb loops and hand pockets make it easier to secure straps during use.



Lightweight

Designed to make all-day use more easy and enjoyable with a soft inner lining. The thin, breathable belt material conforms to user's body while preventing any build-up of moisture.



Easily Removable

Included velcro tabs make it easy to remove brace when finished using.

Distributed by

**vive
health**

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

(239) 220-5772 [vivehealth.com](http://www.vivehealth.com) service@vivehealth.com

8955 Fontana Del Sol Way, Naples, FL 34109, USA